

# RNA-seq


## Quality Control

Aula 02

<https://ttdorres.github.io/transcriptomics/>


# RNA-seq workflow

## Quality Control



# RNA-seq workflow

## Dataset



# RNA-seq workflow

## Dataset

- Calb Female & Larvae
- Cbez Female & Larvae
- Chom Female & Larvae
- Cmac Female & Larvae
- Cmeg Female & Larvae
- Lexi Female & Larvae

# RNA-Seq: initial processing

## Quality check

1. Go to home folder and create a new folder called `rnaseq`

```
cd ~  
mkdir rnaseq
```

2. Go to the `rnaseq` folder and create a new folder called `01-RawReads`

```
cd rnaseq  
mkdir 01-RawReads
```

3. Go to the `01-RawReads` folder

```
cd rnaseq  
mkdir 01-RawReads
```

4. Download the files (R1 and R2) for your sample from the link indicated on the course page ("save link as")

# RNA-Seq: initial processing

## Quality check

5. Check the md5sum of each file

6. View the first 10 lines of each file

# RNA-Seq: initial processing

## Quality check

5. Check the md5sum of each file

```
md5sum Chom_R1.fastq.gz  
md5sum Chom_R2.fastq.gz
```

6. View the first 10 lines of each file

# RNA-Seq: initial processing

## Quality check

5. Check the md5sum of each file

```
md5sum Chom_R1.fastq.gz  
md5sum Chom_R2.fastq.gz
```

6. View the first 10 lines of each file

```
head Chom_R1.fastq.gz  
head Chom_R2.fastq.gz
```



# RNA-Seq analyses

## Sequence formats

### Fasta



The diagram illustrates the Fasta format structure. It shows three entries, each consisting of a header line starting with a greater-than sign (>) and a sequence line. The labels 'Header' and 'Sequence' are placed to the left of each entry, with lines and dots connecting them to the corresponding parts of the Fasta entries.

```
>VIT_201s0011g03530.1
AATTAAGCATAAATACTCACTCTTACCCCCTTATTTTCTTATCTCTCATCACTTTTGGTGCGAAG
GACCATGAGAACAAGCTGCAATGGGTGTAGGGTTCTTCGCAAGGCATGCAGCCAAGACTGCATCA


>VIT_201s0011g03540.1
CAGGTAGCGTGAAGTTAAACCCTAGCGCTTTAGACAAACAGCTGTAGTCACCGCCCACAAACACC
AGCCTCTGAGACACCACCTCAAACCTTTCCACTTAAATACACATCCCTCACACCCTTTTCAATTC

>VIT_201s0011g03550.1
CATGCAAAGCTGAACGCGATGCTGTGATTGGTGGTAAGTGGTAGTTGAGTAAATTTGACAGTGAA
GCCGAAATGGTAAAAGACTAAGGCTAGAAGTAGAATACCACTGTTCTTCTCATCACGTGGGCCCA
```

# RNA-Seq analyses

## Sequence formats

### Fastq



# RNA-Seq analyses

## Usual steps

- Initial Processing
- Assembly (model vs non-model)
- Annotation of transcripts
- Alignment/pseudoalignment reads to a reference
- Differential gene expression

# RNA-Seq analyses

## Initial Processing

- Demultiplexing
- Removing adapters
- Trimming
- Kmer Normalization

# RNA-Seq: initial processing

## Demultiplexing

Multiplexing is a common strategy due to the sheer amount of data a sequencer is able to generate.

After sequencing, each read may be traced back to its original sample using the index sequence and binned accordingly.

- FASTX-Toolkit
- Stacks
- biopieces

# RNA-Seq: initial processing

## Removing adapters

If cDNA insert sizes are sufficiently small and sequencing read lengths sufficiently long, it is possible to generate sequencing reads that contain a portion of adapter sequence at the 3'-end.

- FASTX-Toolkit
- Stacks
- biopieces
- trimmomatic

# RNA-Seq: initial processing

## Trimming

Reads likely to contain multiple sequencing errors provide less biological information and are expected to hinder assembly and alignment.

Some analyses are more tolerant of error than others. For example, de novo assembly requires much cleaner reads than alignment to a reference genome.


- FASTX-Toolkit
- Stacks
- biopieces
- trimmomatic

**How do we know our sequences need trimming?**




# RNA-Seq: initial processing

## Quality check



# RNA-Seq: initial processing

## Conda environments



# RNA-Seq: initial processing

## Quality check

1. In the folder with the sequences, activate `rnaseq` environment

```
conda activate rnaseq
```

2. Call fastqc

```
fastqc
```

# RNA-Seq: initial processing

## Quality check: FASTQC

The main functions of **FASTQC** are:

- Import of data from BAM, SAM or FastQ files
- Providing a quick overview to tell you in which areas there may be problems
- Summary graphs and tables to quickly assess your data
- Export of results to an HTML based permanent report
- Offline operation to allow automated generation of reports without running the interactive application

# RNA-Seq: initial processing

## Running FASTQC

1. GUI: Call fastqc

```
fastqc
```

2. Individual sequence

```
fastqc SpeciesF_R1.fastq.gz
```

3. Bulk (using wildcard `*`)

```
fastqc *.gz
```

4. Using multiple threads

```
fastqc *.gz -t 8
```

# RNA-Seq: initial processing

## Running FASTQC

Using multiple threads

```
fastqc *.gz -t 8
```

### In my computer:

- it took ~5s for two files:
  - ChomF\_R1.fastq.gz with 6.1MB
  - ChomF\_R2.fastq.gz with 6.0MB
- for 2 files with ~2GB each (20million reads), it took ~5min

# RNA-Seq: initial processing

## Running FASTQC

### In your computer:

Use fastqc in a single file:

```
fastqc Species_R1.fastq.gz
```

# RNA-Seq: initial processing

## Running FASTQC

### In your computer:

Use fastqc in a single file:

```
fastqc Species_R1.fastq.gz
```

FASTQ will generate two files:

- Species\_R1\_fastqc.html
- Species\_R1\_fastqc.zip



# RNA-Seq: initial processing

## FASTQC report

- Open `Species_R1_fastqc.html` (double-click):


# RNA-Seq: initial processing

## FASTQC report

- Basic Statistics
- Per base sequence quality
- Per tile sequence quality
- Per sequence quality scores
- Per base sequence content
- Per sequence GC content
- Per base N content
- Sequence Length Distribution
- Sequence Duplication Levels
- Overrepresented sequences
- Adapter Content

# RNA-Seq: initial processing



## Per base sequence quality



# RNA-Seq: initial processing

## Per base sequence quality



### Overclustering



# RNA-Seq: initial processing

## Per base sequence quality

### Issues with the sequencing instruments



# RNA-Seq: initial processing


## Per base sequence content

### RNA-seq

Always gives a FAIL for RNA-seq data. This is because the first 10-12 bases result from the 'random' hexamer priming that occurs during RNA-seq library preparation. This priming is not as random as we might hope giving an enrichment in particular bases for these initial nucleotides.


# RNA-Seq: initial processing

## Per sequence GC content



# RNA-Seq: initial processing

## Sequence Duplication Levels





# RNA-Seq: initial processing

## Trimming and filtering: Trimmomatic

- Cut adapter / other illumina-specific sequences from the reads;
- Cut bases off the start of a read, if below a threshold quality
- Cut bases off the end of a read, if below a threshold quality
- Drop the read if it is too short (below 25 bases)
- Drop low quality reads

# RNA-Seq: initial processing

## Trimming and filtering: Trimmomatic

1. In `rnaseq`, create folder `02-FilteredReads`

```
mkdir 02-FilteredSeqs  
cd 02-FilteredReads
```

2. Save the file `TruSeq3-PE.fa` within `02-FilteredReads`
3. Run trimmomatic in the 'rnaseq' environment

```
trimmomatic  
Usage: PE [-version] [-threads <threads>] [-phred33 |  
-phred64] [-trimlog <trimLogFile>] [-summary  
<statsSummaryFile>] [-quiet] [-validatePairs]  
[-basein <inputBase> | <inputFile1> <inputFile2>]  
[-baseout <outputBase> | <outputFile1P>  
<outputFile1U> <outputFile2P> <outputFile2U>]  
<trimmer1>...
```

# RNA-Seq: initial processing

## Running Trimmomatic

Option	Function
SE or PE	Reads are single end or paired end
ILLUMINACLIP	Perform adapter removal
SLIDINGWINDOW	Perform sliding window trimming
LEADING	Cut bases off the start of a read, if below a threshold
TRAILING	Cut bases off the end of a read, if below a threshold
CROP	Cut the read to a specified length.
HEADCROP	Cut a specified number of bases from the start of the read
MINLEN	Drop an entire read if it is below a specified length
TOPHRED33	Convert quality scores to Phred-33
TOPHRED64	Convert quality scores to Phred-64

# RNA-Seq: initial processing

## Running Trimmomatic

3. Copy the command to run trimmomatic
4. Change the name of the `SPECIES` variable

```
trimmomatic PE -threads 8 -phred33 \  
~/rnaseq/01-RawReads/Chom_R1.fastq.gz \  
~/rnaseq/01-RawReads/Chom_R2.fastq.gz \  
Chom_R1_paired.fastq Chom_R1_unpaired.fastq \  
Chom_R2_paired.fastq Chom_R2_unpaired.fastq \  
ILLUMINACLIP:../02-FilteredReads/TruSeq3-PE.fa:2:30:10:2:\  
SLIDINGWINDOW:4:20 \  
LEADING:3 TRAILING:3 MINLEN:36 >Chom_R1-report.txt
```

# RNA-Seq: initial processing

## Running Trimmomatic

3. Copy the command to run trimmomatic
4. Change the name of the `SPECIES` variable

```
trimmomatic PE \  
    -threads 8 \  
    -phred33 \  
    -summary ${SPECIES}-report.txt \  
    ~/rnaseq/01-RawReads/${SPECIES}_R1.fastq.gz \  
    ~/rnaseq/01-RawReads/${SPECIES}_R2.fastq.gz \  
    ${SPECIES}_R1_paired.fastq ${SPECIES}_R1_unpaired.fastq \  
    ${SPECIES}_R2_paired.fastq ${SPECIES}_R2_unpaired.fastq \  
    ILLUMINACLIP:../TruSeq3-PE.fa:2:30:10:2:True \  
    SLIDINGWINDOW:4:20 \  
    LEADING:3 \  
    TRAILING:3 \  
    MINLEN:50
```

# RNA-Seq: initial processing

## Running Trimmomatic

Input Read Pairs: 100000

Both Surviving: 91421 (91.42%)

Forward Only Surviving: 4513 (4.51%)

Reverse Only Surviving: 2034 (2.03%)

Dropped: 2032 (2.03%)

TrimmomaticPE: Completed successful

# RNA-Seq: initial processing

## Running FASTQC on trimmed data

### In your computer:

Use fastqc in a single file:

```
fastqc Species_R1_paired.fastq
```



FASTQ will generate two files:

- Species\_R1\_paired\_fastqc.html
- Species\_R1\_paired\_fastqc.zip

# RNA-Seq: initial processing

## Running FASTQC on trimmed data

### Compare before and after reports





# RNA-seq workflow

## Quality Control

